

Water Bath  
Slide Dryer

AMOS<sup>®</sup>  
SCIENTIFIC

BETTER HISTOLOGY SOLUTION  
BETTER LIFE FROM HEALTH



## AWB 210 WATER BATH

### Technical data of AWB 210 Water Bath

- Water trough dimensions: 230 × 180 × 56 mm
- Temperature control accuracy: ± 1°C
- Temperature range: ambient temperature ~75°C
- Rated voltage: 220V ± 10% 50/60 Hz or 110V ± 10% 50/60 Hz
- Nominal power: 350W
- Outside dimensions (W × D × H): 328 × 300 × 85 mm
- Weight: about 2.4kg
- Capacity : 2.3 L

