



BETTER HISTOLOGY SOLUTION
BETTER LIFE FROM HEALTH



ASD 220 SLIDE DRYER

Technical Specifications

- Temperature range: ambient temperature ~90°C
- Temperature control accuracy: $\pm 1^\circ\text{C}$
- Capacity : can load 60 slides simultaneously
- Rated voltage:
220V $\pm 10\%$, 50/60 Hz / 110V $\pm 10\%$, 50/60 Hz
- Nominal power: 330W
- Outside dimensions(mm): 330(W) \times 305(D) \times 87(H)
- Weight: 2.75kg

